

HYPER-a

HYPER-SPECTRAL ABSORPTION INSTRUMENT

First *In Situ* Measurements with Hyper-a

- Measurements with Sequoia Scientific's Hyper-a and Sea-Bird Scientific's ac-s collected in East Sound, WA on May 20, 2024, during a period of strong *Noctiluca scintillans* bioluminescence.
- The comparison is not strictly one-to-one, as the two instruments sampled slightly different water volumes.
- The ac-s remain the only mature, routinely available commercial in situ hyperspectral absorption meter.
- Sequoia's Hyper-a, commercially available since 2024 provides much higher spectral resolution, revealing narrow spectral features missed by the ac-s. And extends into the UV (300 nm – 710 nm).





Figure 1. Comparison of hyperspectral absorption measurements (a_{pg} , left y-axis) and their second derivatives (right y-axis) collected on May 20, 2024 in East Sound, Orcas Island, WA, USA. Data from two instruments: Hyper-a (red line) and ac-s (blue line with markers). Prominent absorption features are seen at 379, 411, 439, 462, 493, 589, 618, 637, and 671 nm, corresponding to known pigment signatures in *Chaetoceros socialis* and *Pyrocystis fusiformis* populations present during the *Noctiluca scintillans* bloom. Second-derivative spectra are overlaid to highlight band minima.



Both instruments capture the main absorption envelope, confirming consistent total pigment absorption and show overall shape agreement. Hyper-a has a spectral resolution advantage: Its dual-spectrometer, hyperspectral system (~2 nm resolution) detects narrow pigment features that the ac-s (with ~4–10 nm spacing) does not resolve.

- The 440–460 nm maximum is diagnostic of peridinin-containing dinoflagellates, confirming Noctiluca's carotenoid signature. The broad peridinin band (440–470 nm) appears smoother and slightly stronger in Hyper-a, while ac-s underestimates the peak amplitude.
- Hyper-a resolves finer structure: Around 500–550 nm, Hyper-a shows subtle curvature and a defined minimum, while ac-s produces a flatter, featureless response.
- The valley around 550 nm represents the green window where peridinin and chlorophyll absorb weakly — this gives Noctiluca its orange-brown appearance.
- The 675 nm chlorophyll-a peak confirms the presence of photosynthetic pigments (though Noctiluca is mixotrophic, not a strong chlorophyll producer).

Data, figures, and descriptions courtesy of Alberto Tonizzo, Chris Strait, Mike Twardowski, Sunstone Scientific LLC, and Florida Atlantic University – HBOI, USA. The Hyper-a was developed with NASA SBIR funds on contracts #80NSSC20C0472 and 80NSSC21C0491 to Sequoia Scientific, Inc. with Sunstone Scientific, LLC as subcontractor.

More information: info@sequoiasci.com • <https://www.sequoiasci.com/product/hyper-a/>




Figure 2. Email us for more Hyper-a information!




Figure 3. Hyper-a webpage link

Sequoia Scientific, Inc.
2700 Richards Road, Suite 107
Bellevue, WA 98005 USA
Tel +1 (855) 753-3313
email info@SequoiaSci.com
www.SequoiaSci.com



FLORIDA ATLANTIC UNIVERSITY
Harbor Branch
Oceanographic Institute